

## Gastric Parietal Cell IgG Serum

Gastric Parietal Cell IgG Serum is offered at the Immunopathology laboratory at the main campus.

The new test is a FDA-approved semi-quantitative assay to be used as an aid in diagnosing autoimmune gastritis and pernicious anemia. The numerical values (ELISA units) are only reported for informational purposes and may not correlate with presence of clinically overt disease or disease severity. The new assay targets purified gastric parietal cell H<sup>+</sup>/K<sup>+</sup> ATPase which is the main target in autoimmune gastritis.

This test offers improved sensitivity and specificity compared with the IFA. Gastric parietal cell IgG may be present years before clinical manifestations and complications such as pernicious anemia. The seroprevalence increases with age and has been reported in 2% of normal population younger than 20 years and in 16% of population older than 60 years; this is particularly the case in those with concomitant Hashimoto thyroiditis, Sjögren's syndrome, and concomitant *Helicobacter pylori* infection, among other conditions, therefore, clinical correlation is highly recommended.

### Test Overview

<b>Test Name</b>	Gastric Parietal Cell IgG Serum
<b>Ordering Mnemonic</b>	PARIES
<b>Methodology</b>	Enzyme-linked immunosorbent assay (ELISA) (aka EIA)
<b>Specimen Requirements</b>	Primary container SST (Gold) 1.0 mL; Minimum 0.3 mL
<b>Stability</b>	After separation from cells: 24 hours (Ambient), 7 days (Refrigerated), 14 days (frozen ≤ -20°C)
<b>Clinical Information</b>	As an aid in diagnosing autoimmune gastritis
<b>Reference Range</b>	Negative; ≤20
<b>CPT Code</b>	83516
<b>Days Performed</b>	Tuesday, Friday

### Technical Information Contact:

Matthew Dee, C(ASCP)  
Supervisor, Immunopathology  
216.445.0511

### Clinical Information Contact:

Kamran Kadkhoda, PhD, D (ABMM), D(ABMLI)  
Medical Director, Immunopathology Laboratory  
kadkhok@ccf.org